

CaviWipesXL™

extra large
disinfecting
towelette

- TUBERCULOCIDAL
- BACTERICIDAL
- VIRUCIDAL*
- FUNGICIDAL

Metrex

Net Contents: 1 Pre-Saturated Towelette (12"x10")

Packed For: METREX RESEARCH CORPORATION®
28210 Wick Road, Romulus, MI 48174

CaviWipesXL Disinfecting Towelette is a non-woven disposable towelette pre-saturated with CaviCide®. **CaviWipesXL** is intended for use in health care settings such as hospitals, laboratories, clinics, dental offices, ophthalmic offices and veterinary facilities. **CaviWipesXL** is effective against the following microorganisms on hard, non-porous surfaces when used as directed:

- *Mycobacterium tuberculosis* var: *bovis* (BCG)
- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Salmonella choleraesuis*
- *Trichophyton mentagrophytes*
- Methicillin Resistant *Staphylococcus aureus* (MRSA)
- Vancomycin Resistant *Enterococcus faecalis* (VRE)
- Hepatitis B Virus (HBV)*
- Herpes simplex virus types 1 and 2*
- Human Immunodeficiency Virus (HIV-1)*

ACTIVE INGREDIENTS:

Diisobutylphenoxyethoxyethyl dimethyl benzyl ammonium chloride0.233%
Isopropanol.....14.300%

INERT INGREDIENTS:.....85.467%
TOTAL100.000%

PRECAUTIONARY STATEMENTS: Hazards to Humans and Domestic Animals CAUTION

Harmful if absorbed through the skin. Causes moderate eye irritation. Avoid contact with eyes, skin, or clothing. Wash hands before eating, drinking, chewing gum, using tobacco, or using the toilet. Remove contaminated clothing and wash clothing before reuse.

EPA Reg. No. 46781-8 • EPA Est. No. 71979-SC-001 • 56952-WI-001

KEEP OUT OF
REACH OF
CHILDREN
CAUTION

See back panel for
additional precautionary
statements.