

# ORACLE

# 1

## ONE STEP•ONE MINUTE CLEAN AND DISINFECT

Virucidal\* • Bactericidal • Tuberculocidal • Fungicidal

### ACTIVE INGREDIENTS

Ethyl Alcohol .....	72.50%
Didecyl Dimethyl Ammonium Chloride .....	0.33%
OTHER INGREDIENTS: .....	27.17%
Total .....	100.00%

## KEEP OUT OF REACH OF CHILDREN CAUTION

See side panel for additional precautionary statements

### PRECAUTIONARY STATEMENTS

Hazards to Humans and Domestic Animals:

**CAUTION:** Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the toilet.

**FIRST AID: IF IN EYES:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice.

Have the product container or label with you when calling the poison control center or doctor, or going for treatment. For emergency information concerning this product, call the National Pesticides Information Center (NPIC) at 1-800-858-7378 seven days a week, 6:30 am to 4:30 pm Pacific Time [NPIC Web site: [www.npic.orst.edu](http://www.npic.orst.edu)]

**PHYSICAL OR CHEMICAL HAZARDS:** Flammable. Keep away from heat or open flame.

**Oracle is intended for use on hard non-porous surfaces in:**

• Hospitals and healthcare settings, including: Ambulances, Anesthesia, CAT Labs, CCU, Central Supply, Clinics, Dental Offices, Dialysis Clinics, Doctor's Office, Emergency Rooms, Home Health

Care, Hospices, Laboratory, Laundry Rooms, Long Term Care Facilities, Newborn Nursery, Nursing Homes, Offices, Operating rooms, Pediatrics, Physical Therapy, Physician's Office[s], Ophthalmic Offices, Radiology, Recovery Room, Rehabilitation, Respiratory Therapy, Ultrasound, Transport Vehicles, Waiting Areas, X-Ray.

- Critical Care Areas: [CCU] Cardiac Care Units, [E.R.] Emergency Rooms, [ICU] Intensive Care Units, [NICU] Neonatal Intensive Care Units, [PICU] Pediatric Intensive care Units, Surgery, Surgical Prep, and Surgical Recovery Units.
- Childcare Facilities, Daycare Centers, Nurseries.
- Households, Homes, Residences and/or Bathrooms.
- Institutional Bathrooms.
- Animal Housing/Holding Facilities, Animal Hospitals, Animal Life Science Laboratories, Animal Transportation Vehicles, Breeding Establishments, Grooming Establishments, Kennels, Pet Shops, Tack Shops, Veterinary Clinics, Veterinary Examination Rooms, Veterinary Operating Rooms, Veterinary Waiting Rooms, Veterinary X-Ray Rooms, Zoos.

Net Contents:

32 Fl. Oz.



Sold By: Share Corp.

P.O. Box 245013

Milwaukee, WI 53224

(800) 776-7192 • [www.sharecorp.com](http://www.sharecorp.com)

### DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body.

**CLEANING:** Wet surface sufficiently and clean with a mop sponge, cloth, nonwoven wipe or other suitable cleaning instrument.

**DISINFECTING:** To disinfect hard, non-porous surfaces, thoroughly wet the surface to be treated. Treated surface must remain visibly wet for one minute to achieve complete disinfection of all pathogens listed on this label. For visibly soiled surfaces, clean before following disinfecting instructions. Allow surfaces to air dry.

**\*VIRUCIDAL:** Oracle is an effective virucide on hard, non-porous surfaces, in the presence of an organic load, against *Adenovirus type 5* [ATCC VR-5], *Avian Influenza A Virus* [ATCC VR-798], *Cytomegalovirus* [ATCC VR-538], *Enterovirus EV-D68* [ATCC VR-1826], *Hepatitis B Virus*, *Hepatitis C Virus*, *Herpes Simplex Virus Type 1* [ATCC VR-539], *Human Coronavirus* [ATCC VR-1558], *Human Immunodeficiency Virus* [HIV-1] [AIDS], *Influenza Virus type A2* [ATCC VR-544], *Influenza B* [ATCC VR-102], *Norovirus*, *Poliovirus serotype 1* [PV-1] [ATCC VR-1562], *Respiratory Syncytial Virus* [ATCC VR-26], *Rhinovirus* [ATCC VR-340], *Rotavirus Strain WA* [ATCC VR-2018], [and] *Vaccinia virus* [ATCC VR-119], when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

**KILLS HIV-1, HBV and HCV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings [Hospitals, Nursing Homes] or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 [HIV-1] [associated with AIDS], Hepatitis B Virus [HBV], and Hepatitis C Virus [HCV].

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV, and HCV FOR SURFACES OR OBJECTS SOILED WITH BLOOD OR BODY FLUIDS.**

**PERSONAL PROTECTION:** When handling items soiled with blood or body fluids, use disposable latex gloves, gowns, masks and eye coverings.

**CLEANING PROCEDURES:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of Oracle.

**CONTACT TIME:** Allow surface to remain wet for 1 full minute.

**DISPOSAL OF INFECTIOUS MATERIALS:** Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

**FUNGICIDAL:** Oracle is an effective fungicide against *Trichophyton mentagrophytes* [ATCC 9533], causative agent of Athlete's Foot Fungus, when used as directed in bathrooms, shower stalls, locker rooms, exercise facilities or other surfaces commonly contacted by bare feet, and *Candida albicans* [ATCC 10231]. The surface must remain visibly wet for 1 minute. Allow surfaces to air dry.

**BACTERIA:** Oracle is an effective disinfectant on hard non-porous surfaces, in the presence of an organic load, against *Escherichia coli* [E. coli] 0157:H7 [ATCC 35150], *Pseudomonas aeruginosa* [ATCC 15442], *Salmonella enterica* [Salmonella] [ATCC 10708], *Staphylococcus aureus* [ATCC 6538], *Streptococcus pyogenes* [Strep] [ATCC 19615], *Klebsiella pneumonia* [ATCC 4352], *Listeria monocytogenes* [ATCC 984], *Campylobacter jejuni* [ATCC 29428], and *Burkholderia cepacia* [ATCC 25416], when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

### DRUG RESISTANT AND MULTI-DRUG RESISTANT BACTERIA

Oracle is an effective disinfectant on hard non-porous surfaces, in the presence of an organic load, against *Acinetobacter baumannii* – Multi-Drug resistant, *Carbapenem resistant Enterobacter cloacae* (CRE) (KPC) [Multi-Drug resistant], *Carbapenem resistant Enterobacter cloacae* (CRE) (New Delhi) [Multi-Drug resistant], *Enterococcus faecium* – Vancomycin resistant (VRE), *Escherichia coli* – ESBL Producing [Multi-Drug resistant], *Klebsiella pneumoniae* – ESBL Producing [Multi-Drug resistant], *Staphylococcus aureus* – Methicillin resistant (MRSA), *Staphylococcus aureus* – Vancomycin resistant (VRSA), and *Staphylococcus aureus* Panel – Vancomycin intermediate (VISA) , when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

**TUBERCULOCIDAL:** Oracle is an effective disinfectant on hard non-porous surfaces, in the presence of an organic load, against *Mycobacterium bovis* [ATCC 35743], when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

### ANIMAL PATHOGENS

Oracle is an effective disinfectant on hard non-porous surfaces against *Avian Influenza A Virus*, *Bordetella bronchiseptica*, *Bovine Diarrhea Virus* (BVDV), *Canine distemper Virus*, *Canine parvovirus*, *Duck Hepatitis B Virus* (BVDV), *Feline calicivirus*, and *Feline leukemia Virus*.

**Veterinary Clinics, Animal Life Science Laboratories, Pet Shops:** Remove all animals and feed from premises, vehicles, and enclosures. Remove all litter and manure from floors, walls, and surfaces of facilities and fixtures occupied or traversed by animals. Empty all troughs, racks, and other feeding and watering appliances. Thoroughly clean all surfaces with soap or detergent and rinse with water. Wet surfaces for 1 minute. Immerse all handling restraining equipment as well as forks, shovels, and scrapers used for removing litter and manure. Ventilate buildings, vehicles, and other closed spaces. Do not house animals or employ equipment until treatment has been absorbed, set, or dried. Thoroughly scrub treated feed troughs, racks and other feeding and watering appliances with soap or detergent, then rinse with potable water before reuse.

**EFFICACY TESTS HAVE DEMONSTRATED THAT ORACLE IS AN EFFECTIVE, VIRUCIDE, BACTERICIDE, TUBERCULOCIDE AND FUNGICIDE IN THE PRESENCE OF 5% ORGANIC SOIL.**

Oracle will not harm most surfaces: including stainless steel, aluminum, other hard and soft metals, sealed quartz or other solid man-made counter top materials, most plastics, fiberglass, vinyl, nylon, and other non-porous upholstery materials used in a medical, veterinary or dental environment, and most painted surfaces. For plastic and painted surfaces, spot test on an inconspicuous area before use. Not recommended for use on acrylic plastics or coated optical surfaces. Do not use on unpainted wood, or leather. Not for use on people or pets.

## STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal.

**PESTICIDE STORAGE:** Always store this product and all pesticides in the original container. If a leaky container must be contained within another, mark the outer container to identify the contents. Protect from extreme heat or cold.

**PESTICIDE DISPOSAL:** Waste resulting from the use of this product must be disposed of on site or at an approved waste disposal facility.

**CONTAINER DISPOSAL:** Nonrefillable container. Do not reuse or refill this container. Offer for recycling if available, or discard empty container in trash.

EPA Reg. No. 88494-3-11547 EPA Est. No. 10088-WI-2

MADE IN USA